

Analytical Balance

220g x 0.1mg (Economy)

MODEL: MA-220X



SALIENT FEATURES:

- Strong anti-shaking monolithic electromagnetic sensor
- 18 standard units, with character mark
g, kg, ct, lb, oz, ozt, tlh, tls, GN, dwt, mg, /lb, tlc
mom, tol, bat, MS
- In-built density weighing function
- Dynamic animal weighing function

SPECIFICATIONS:

Model	MA-220X
Capacity	220g
Make	Galaxy
Readability	0.0001g (0.1mg)
Linearity	± 0.0003g
Minimum Load	0.01g
Display & Pan Size	LCD, Pan: Ø 90 mm
Response Time	4-6 sec
Calibration	External Calibration
Tare Range	Full
Warm Up time	30 min
Units	g, ct, ozt, oz, GN, dwt, lb & 11 more
Weighing Chamber Height	200 mm
Class	I
Interface	RS-232C
Power Supply	AC 220V / 50 Hz
Dimensions & Weight	200 x 290 x 310 mm & 9.5 kg

STANDARD

ACCESSORIES:

- Draft Ring
- Pan
- Power Cord
- Instruction Manual