


MODEL: HE-MAXI



SPECIFICATIONS:

Model	HE-MAXI
Tray size (cm)	16 x 20, 16 x 15
Sample Combs	1 x 17 wells (1mm)
	1 x 34 wells (1mm)
	1 x 17 wells (1.5mm)
	1 x 34 wells (1.5mm)
Comb thickness	1.0,1.5 mm
Buffer Volume	1000 ml
Qty. of gel made	1 Piece
Size	370 x 210 x 130 (mm)
Weight	2 kgs

STANDARD

ACCESSORIES:

- Main Buffer Tank
- Lid with Electrode
- Different sizes of combs
- Casting Trays