

Mini Centrifuge With timer & speed display

MODEL: MCX-12



SALIENT FEATURES:

- Brushless motor for long life.
- 2 adapters included to use PCR & MCT tubes.
- Automatic motor shut down on opening the cover.
- Compact design.
- LCD display and microprocessor control.

SPECIFICATIONS:

Model	MCX-12
Make	Galaxy
Speed Range	100 to 12000 RPM
Speed Adjustment	Digital
Speed Increment	10 RPM
Speed Control Accuracy	±3%
Relative centrifugal force	7245 g
Timer	1 to 99 minutes
Motor Type	Brushless DC Motor
Motor Power	48 W
Display	LCD
Noise	51-53 dBA
Rotor	MCT 8*2mL, PCR 16*0.2mL
Tubes which can be used	2/1.5/0.5/0.2mL*8, 2mL*16, 0.2mL*16
Adapter	0.5ml, 0.2ml Adapter
Dimensions & Weight	170 x 146 x 119 mm & 1 kg

STANDARD

ACCESSORIES:

- 8 Hole Rotor
- 16 Hole Rotor
- 0.5ml Adapter
- 0.2ml Adapter
- Power Adapter
- Instruction Manual